


# FTP Access to Data

The NMR data can be accessed from a personal or lab computer via a FTP connection. The FTP connection can only be accessed while on the UTSA network.

Filezilla is the recommended FTP client. <https://filezilla-project.org/>



## IP Addresses

DD2 500 – 10.4.5.39

Username – your account username

Password – your account password

Bruker 500 and 300 – 10.4.2.22 port # 54806

Username – UTSAchemistry

Password – NMR\_Bruker

## File Transfer

To transfer the NMR files, drag the files you want from the NMR computer side (right side) to where you want to save them on your computer (left side).